

Mouvlight Wheelchair



- 40 cm (16 inches) - 876001.40

- 42.5 cm - 876001.42

- 45 cm (18 inches) - 876001.45

- 47.5 cm - 876001.48

(On restock - Available from 23/08/2023)

- 50 cm - 876001.50

(On restock - Available from 23/08/2023)

- **Sturdy foldable double cross structure**
- **Third party brakes**
- **Adjustment capacity**
- **Video available below**

Our MouvLight manual wheelchair is a sturdy and versatile wheelchair for hire or purchase. It comes standard with third party brakes, an ergonomic polypropylene handrail, anti-tilt castors and soft PU armrests with wide clothing protection sides. Quick release 24 inch wheels. Equipped with Velcro fastening between backrest and seat. The Mouvlight is a manual folding wheelchair.

Height adjustment of the 24 inch wheels in 3 positions (-2,5 and +2,5 cm). Offset adjustment of the 24 inch wheels in 4 positions (-5 cm, -2.5 cm and +5 cm). Front fork 3 heights (-2.5 and +2.5 cm). Footrest adjustable in 5 positions from 35 to 45 cm from the seat.

additional information: solid wheels

5 seat widths available: 40 / 42.5 / 45 / 47.5 / 50 cm.

Overall dimensions: width 62 / 64,5 / 67 / 69,5 / 72 cm x depth 78 cm without footrest, 102 cm with x height 90 cm. Foldable up to 36 cm wide.

Seat dimensions: width 40 / 42,5 / 45 / 47,5 / 50 cm x depth 42 x height. 47 cm.

Colour: metal grey.

Weight: 15 to 16.5 kg depending on width, fully equipped.

Max weight : 125KG

Variantes :

- 40 cm (16 inches)
- 42.5 cm
- 45 cm (18 inches)
- 47.5 cm
- 50 cm