

Robust wheelchair



- 874002
(On restock - Available from 23/08/2023)

- **Robust and foldable wheelchair for a good price**
- **Seat length 45 cm**

This robust wheelchair is designed for persons with reduced mobility living in medicalised structures.

It is equipped with pull-down armrests, removeable armrests and a folding steel frame.

Rear wheels of 24 inches and front wheels of 20 cm.

Seat width 45 cm.

aLuminium colour

Over-all dimensions: width 63 x depth 103 x height 91 cm.

Seat dimensions: width 45 x depth 40 x height 50 cm.

Backrest height 44 cm

Weight 18.5 kg

Max weight allowed 110 kg