

Modulolight transfer chair



- 42.5 cm - Metal Grey - 874006.G
- 42.5 cm - Metal Red - 874006.R
- 45 cm (18 inches) - Metal Grey - 874007.G
- 45 cm (18 inches) - Metal Red - 874007.R
- 47.5 cm - Metal Grey - 874001.G
- 47.5 cm - Metal Red - 874001.R

- **Standard self-locking third-party brakes**
- **Double cross-brace as standard**
- **Flip-up armrests and clothes guard as standard**

Our Modulolight transfer chair is ideal for daily transfers at home or in medicalised structures. It has a very sturdy aluminium frame with double crossbars that can accommodate people weighing up to 125 kg. Its large soft wheels (30 cm at the back and 20 cm at the front) offer a comfortable outdoor walks. Self-locking third party brakes as standard. Foldable armrests.

3 widths available: 42.5 cm, 45 cm and 47.5 cm.

Overall dimensions: width. 58.5 / 61 / 63.5 x depth. 97 x height. 92 cm.

Seat dimensions: width 42.5 / 45 / 47.5 x depth 41 x height. 49 cm. Backrest height 43 cm.

Folded dimensions (including folded breakable backrest): width 27 x depth 65 x height 72 cm.

Weight 12.4 kg.

Resistance 125 kg.

2 colours available: metal grey and metal red

Variantes :

- 42.5 cm
- 45 cm (18 inches)
- 47.5 cm

- Metal Grey
- Metal Red