

Tiger Light Rollator



- 826269

- **Handy and very light**
- **Height-adjustable backrest**
- **No need to remove the bag to fold it up**

Our Tiger Light walker allows you to walk and shop with maximum support and safety. It has a very practical cane holder, a useful walkway for the city and a height-adjustable comfort backrest. Its seat is upholstered for extra seating comfort. Its front wheels are removable for more compact storage. Front wheels 20 cm, rear wheels 17 cm.

Overall dimensions: width 64 x depth 73 x height 86 / 98.5 cm in 6 positions. Weight 8 kg.

Seat dimensions: width 45 x depth 25 x height 54 cm.

Maximum supported weight : 130 kg