

PROFILO BATH SEAT

Reference 812138



Description :

The Profilo bath seat lets you wash comfortably and safely while seated. You can choose between 3 seat heights: the high seat (see photo) is particularly suitable for people with reduced mobility, who will be able to straddle the edge while already seated. The lower seats allow you to sit below the level of the ledges. Made to the highest quality standards, it is easy to fit onto the rim of your bathtub, without the need for tools, and holds firmly in place without moving. The seat's side supports are adjustable in width, so the profilo fits most bath models.

Maximum user weight : 110kg

Certification and standards: EC Declaration of Conformity in accordance with Regulation (EU) 2017/745, EN ISO 14971:2019, PN-EN 1041:2010, PN-EN 12182:2012, EN ISO 14001:2015, EN ISO 13485:2016

Cleaning :

Clean with non-abrasive, non-aggressive and non-corrosive household products.
Removable comfort cushion for easy cleaning.

Packaging dimensions :

Box : 510 x 455 x 116 mm


Gross weight : 3,309kg

Net weight : 2,701 kg

Made in Italy.

Distributor : Identités

Exploded view



| Number | Part | Quantity | Material |
|--------|----------------|----------|-----------------|
| 1 | Seat | 1 | PP + 20% talc |
| 2 | Seat support | 2 | PP + 20% talc |
| 3 | Covers | 8 | PP |
| 4 | Joints | 4 | PA6 + 30% FG |
| 5 | Covers | 4 | PP+20% Talc |
| 6 | Extrusion bars | 2 | Aluminium |
| 7 | Screws | 8 | Stainless Steel |
| 8 | Washers | 4 | Stainless Steel |
| 9 | Screws | 4 | Stainless Steel |
| 10 | Anti-slips | 6 | PVS SOFT |
| 11 | Spacers | 4 | PP |